

## 94400 - 94416 ORASOL® Dyestuff

These dyestuffs are soluble in most solvents, are very brilliant and have a good stability.  
All these dyestuffs are easily soluble in a mixture of 80% ethyl alcohol (70800), 5% water and 15% Methoxypropanol PM (70920). For most applications 1-5 g are sufficient for 1 liter.

### Lightfastness and Solubility:

<i>Product name</i>	<i>Lightfastness 5 <math>\mu</math>m</i>	<i>Lightfastness 10 <math>\mu</math>m</i>	<i>Solvent mixture g/l</i>
94400 Yellow 4GN, greenish	4	5	900
94402 Yellow 152, sun gold (old Yellow 2 GLN)	5	7	500
94404 Yellow 190, yellow-orange (old Yellow 2 RLN)	4	5-6	1000
94406 Orange 247 (old Orange G)	4	5	150
94408 Orange 272, orange-red (old Orange RG)	5	6	150
94410 Brown 324 (old Brown 2RL)	4-5	5	150
94412 Red 395, red-violet (old Red BL)	4	5	250
94414 Blue 825, turquoise dark (old Blue GN)	4-5	6	450
94416 Blue 855, turquoise light (old Blue GL)	5	6	25